

Etac Transfer Boards

Mind the gap...

A transfer board is intended to bridge shorter gaps when performing seated transfers. The board's low friction upper surface allows the user to transfer independently or with assistance. Transfer boards are used when transferring between bed and wheelchair / shower chair, between wheelchair / chair / toilet chair or to and from the car.

All Etac Transfer Boards have a non-slip surface on the underside. Combining a transfer board with a glide cushion provides even lower friction and is recommended when the transfer involves bare skin, for instance to and from a shower chair or toilet.

Etac 3B-Board

Etac 3B-Board is a stable and simple transfer board with a traditional curved shape. The board can be placed either curve forward or curve backward depending on the transfer situation and surrounding obstacles, eg: the wheel of the wheelchair. 3B-Board's length allows transfer across wider gaps.




Etac E-Board

Etac E-Board is a thin and flexible transfer board that allows easier application and removal.

E-Board Wings have wings along the sides, which can be bent up or down. This provides more sturdiness and stability, and the upward wing is used to protect against the wheel of the wheelchair. The 45 cm board is more suitable for active and independent users, while the 60 cm boards suit the most common transfer situations.

E-Board (without wings) is smaller and lighter. As it is easily carried, it is suitable for the more active and independent user. The 60 cm boards suit the most common transfer situations while the 75 cm board is suitable for longer transfer, eg. in and out of a car.



Intended Use

- Transfer from bed to wheelchair / shower chair
- Transfer in and out of the car
- Transfer between wheelchair, chair, toilet, shower chair, commode



Description

| | | | | |
|--|------------|-------------|-------------|-------|
| 3B-Board, curved transfer board white | 23 cm 9.1" | 77 cm 30.3" | | IM403 |
| E-Board Small, transfer board black | 20 cm 7.9" | 45 cm 17.7" | 0.5 cm 0.2" | IM410 |
| E-Board Medium, transfer board black | 25 cm 9.8" | 60 cm 23.6" | 0.5 cm 0.2" | IM408 |
| E-Board Medium, transfer board wings black | 33 cm 13" | 60 cm 23.6" | 0.5 cm 0.2" | IM404 |
| E-Board Long, transfer board wings black | 33 cm 13" | 75 cm 29.5" | 0.5 cm 0.2" | IM405 |

Specifications

Polyethylen

max 150 kg max 330 lbs



Width

Length

Height

Item no.

Washing Instructions

Transfer boards can be cleaned with soapy water or disinfectant.

