

EVO CARE CONTACT:

Sales: 07 3355 8000

Fax: 07 3355 5043

Email: sales@evocare.com.au

PO Box 1144, Stafford

Q 4053 Australia

Gel Wheelchair Armrest Pad


Call or email for more information

RH-3 Gel Wheelchair Armboard Pad

The **RH-3 Gel Wheelchair Armboard Pad** is made with high-quality polyurethane gel.

- It provides excellent pressure relief and great comfort.
- The Gel material will not migrate, leak, evaporate, harden, or dry out if the cushion is punctured.
- The Gel is highly resistant to bacteria and to bodily fluids, and can be easily cleaned with warm soapy water or with a disinfectant.
- Dimensions of Gel Pad: 230mm long x 70mm wide x 15mm high
- The armboard pad is manufactured with velcro tape for quick and stable attachment.
- Supplied individually as single armboard pad.

[Click here for a link to the Gel Wheelchair Arm Pads on the Evocare website](#)



 **Evocare**[®]

Our catalogue is our comprehensive and up to date website: <http://www.evocare.com.au/>

Specialist Equipment Sales & Service