

Elbow / Forearm / Canadian Crutches

The forearm cuffs and handgrips are designed to take some of the weight off of the arms during weight bearing.

Learning how to use and adjust forearm crutches properly will ensure that you do not cause yourself further injury or unnecessary muscle strain.



– Double Adjustable

Product Code	Size	Height Adjustable Range (mm)	Hand Grip to Cuff Adj. (mm)	Hand Grip to Tip Adj. (mm)
LM3806L	Large	1080-1390	240-320	840-1070
LM3806M	Medium	950-1290	240-320	710-970
Safe Working Load		100kgs		

- Overall length is measured from top of cuff to tip.
- Hand grip to cuff is measured from top of cuff to top of hand grip.
- Hand grip to tip is measured from top of hand grip to tip.
- Soft touch foam handgrips
- Shaft constructed from strong, lightweight aluminium.
- **ARTG No:** 175494