

Accessories

Clinic | **Care** | Assist



There's nothing like original
Stieglmeyer **ACCESSORIES.**



 **Evocare**
www.evocare.com.au

Ph: 07 3355 8000 Fax: 07 3355 5043
Email: sales@evocare.com.au

**Top quality and modular
OEM ACCESSORIES.**

We not only manufacture beds for the healthcare sector, but also provide complete system solutions that are designed to meet the increasingly specialised needs of healthcare providers.

To this end, we place great emphasis on making modular products that can be optimally adapted to the needs of your facility.

TABLE OF CONTENTS

03	Anti-fall devices
04	Bed extensions
08	Protective covers
10	Accessorised patient lifting pole
12	Lamps
14	Accessories for extra mobility
16	Care accessories
18	Transport and other convenience accessories
20	Additional convenience accessories
22	Bedside table accessories

Safety first: a full range of anti-fall devices

A fall, which only takes a fraction of a second, can have far-reaching health consequences for the resident affected, and can lead to a law suit. Anti-fall devices help to avoid such incidents without unduly restricting residents' freedom of movement.


Our safety mat, which is made of robust artificial leather and foam padding (large picture), abuts the bed rails, and when not in use can be simply folded out of the way and attached to the rails via convenient snaps (figure 1).

Our safety ramp, which features foam and a polyurethane cover, converts what could otherwise be a nasty fall into a rolling movement (figure 2).

Decor: Natural beech R 5311 VV



Small safety mat



Bed extension element



Extra bed length for tall residents.

Our bed extension solution allows standard beds to be lengthened to accommodate unusually tall residents so that they can enjoy healthy and restful sleeps. The product can be readily installed by your handyman.



**BED EXTENSION for nursing beds
The Classiko and Elvido-Rega models.**

This bed extension (figure 1) can be factory-installed, or you can mount the extension at your centre. The extension comprises a bed extension element and rail extension elements (figure 3).


The mattress is extended via a 20 cm-long mattress extension element (figure 2) that is placed on the bed extension element.

When ordering, please indicate the height of the relevant mattress.



Telescoping bed extension

Scope of delivery for the nursing home model	Figure	4 bed rails	2 trim elements	1 bed extension element	1 mattress extension element
Elvido with telescoping bed extension	p. 6				x
Elvido-Brevo with telescoping bed extension	p. 6				x
Elvido-Vano with telescoping bed extension	p. 6				x
Elvido	p. 7, figures 2 and 3	x	x		x
Elvido-Brevo	p. 7, figures 2 and 3	x	x		x
Elvido-Vano	p. 7, figures 2 and 3	x	x		x
Elvido-Rega	p. 4-5	x		x	x
Classiko	p. 4-5	x		x	x
Amara	p. 7, figures 4 and 5	x	x	x	x
Pareno	p. 7, figures 4 and 5	x	x	x	x
Belindo	not shown	x	x	x	x
Arena	not shown	x	x	x	x



Additional bed extension solutions

ELVIDO bed extension

All Elvido models integrate a bed extension element that can be pulled out by simply replacing the bed rails and their trim elements. Installed product: see figure 2. Scope of delivery: see figure 3.

Our Elvido, Elvido-Brevo and Elvido-Vano models allow for delivery of a factory-installed* bed extension with telescoping bed rails and trim elements (large picture). These solutions obviate the need for any additional bed extension elements. The bed can be extended by merely pushing a button. Installation kits are available for all other models.

To increase mattress length, a 20 cm-long mattress element is installed (figure 1). When ordering, please indicate the height of the relevant mattress.

Amara and Pareno bed extensions

The Amara and Pareno bed extensions are composed of a bed extension element and bed rail extensions. Installed product: see figure 4. Scope of delivery: metallic bed extension element and bed rail extension element (figure 5).



Protective covers optimally safeguard residents and equipment.


Our protective covers serve a number of purposes, notably safeguarding residents against injury at the sides, head and foot of the bed.

The covers are made of foam covered with artificial leather. The sides facing the bed have a polyurethane coating and are stitched together. They are also maintenance-friendly, durable, and at the same time agreeably soft and easy to grasp. The covers are elastic and are resistant to abrasion, tearing, and to commercial cleaning agents and disinfectants.

The covers are available in the following formats: (a) **covering for a full-bed rail** consisting of a 200 cm and 220 cm rail (figures 1 and 6), or two 200 cm rails plus a vertical element (figure 2), or three 200 and 220 cm rails (not shown); (b) single covering for a combined rail for segmented rails (large picture); (b) **Mobivit**, an entry-level covering that provides a transitional element between two covers (figure 3); (c) **headboard and footboard covering** with velcro on the back for easy attachment (figure 5).




Rapidly removable bed-rail cushions with a removable polyurethane protective cover and particularly soft and contoured foam (figure 4). Easy to insert on a combined segmented rail or for all-bed rails, with two cushions for complete protection.

Dimensions: 950 x 265 x 100 mm.



Patient lifting pole with grip
and device wire holders,
e.g. for call system





Patient lifting pole – the key to freedom of movement.

Residents need to be able to move around as freely as possible. But some residents need equipment to help them sit up in bed. We have developed various solutions that meet this need.

Our **patient lifting pole** (figure 1) can be swivelled to the right or left side of the bed head. The fixture, which is made of a dual-wall metal pipe, is very stable. Two pins integrated into the end fixture prevent the grip from slipping.

The **grip** is made of robust plastic and features particularly harmonious gripping properties. Its height is readily adjustable via a roller belt (figure 2). The two gripping surfaces feature a comfortable soft-touch surface.

The **wire holders** (figure 3), which clip onto the patient lifting pole, allow for safe and tidy wire installation. They also take the strain off the installed wires and can prevent wire connections and manual operating elements from being yanked out.

A **rope ladder**, which is mounted at the foot of the bed and laid on the bed, enables residents to sit up in bed without using the patient lifting pole (figure 4).

Our plastic **IV bag holder** can be readily attached to the patient lifting pole's end fixture. The holder should be attached adjacent to the front or rear pin of the fixture so as to prevent the holder from rotating (figure 5).



And there was light: our lamp portfolio.

Each resident should have their own personal lamp that does not disturb other residents, as this promotes independence and a sense of well-being.

SOLA

Sola is a versatile halogen reading lamp whose position can be readily adjusted. The lamp is made of fracture-resistant plastic and features a robust, flexible arm. Its ventilation slits and heat dissipation cylinder allow for light that generates a minimal amount of heat (figures 1-3 and 6).

Sola lamp with long arm and illuminated cord switch (figure 1).

Sola lamp with long arm and rocker switch (figure 2).

Sola lamp with short arm (figure 3). The lamp can be mounted on the patient lifting pole via a clamp.

Patient lifting pole with attachment ring to limit the rotation of the Sola lamp (figure 4).

LUMO LED lamp



Maintenance-free reading lamp that provides outstanding, low-energy light without any heat generation. Our Lumo lamp produces bright, focused light ranging up to 1100 lux (large picture).

Bedside table LAMP

Our short-armed Lumo (figure 5) and Sola (figure 6) lamps can be attached to our night table utility bar or the like.



Lumo





Wall spacer

Wall spacer for model	Attached at the headboard	Attached on the side
Elvido	Type 01	Type 01
Elvido-Brevo	Type 01	Type 01
Elvido-Vano	Type 03	Type 01
Elvido-Rega	Type 02	Type 01
Classiko	Type 02	Type 01
Amara	Type 04	Type 01
Amara-Tero	Type 04	Type 01



Additional convenience accessories.

Wall spacers avoid damage to walls and wall sockets. Our wall spacers can be mounted without the use of tools by simply attaching them to the bed frame (large picture).

And speaking of mobility, our **wall protection rollers** avoid wall damage when beds are moved. This applies to both sideways movements (figure 2) and to bed height adjustments thanks to the additionally integrated plastic rollers (figure 3).

If desired, the Elvido bed's entire electrical system can be run using a **24-volt*** external plug-in transformer, which integrates a safety plug that is compatible with all commercial sockets (figure 1). This makes the system completely safe, as the use of 230 volt electricity at residents' beds is avoided. This solution also allows for more rapid adjustment when power consumption is at a low level. The rechargeable battery allows the Elvido bed to be adjusted cordlessly, thus ensuring complete mobility.

The rechargeable battery is available in a set that also includes a charger and holder (figure 4).

Rechargeable-battery holder for nursing beds (figure 5).



Additional convenience accessories.

Our accessories allow essential care elements to be mounted on nursing beds in a stable and practical manner so that residents can be provided with the necessary care. These accessories are compatible from a functionality and design standpoint and offer optimal ergonomics and handling.

With the Elvido bed that integrates supply rails, our **all-purpose drainage-bag holder** can be inserted into the fixture that is provided for this purpose (figure 1).

Our **urine bag holder** can be used to insert a drainage tube between bed rails. The holder is attached to the lower rail, thus protecting the tube (figure 2).

Our urine bag holder (figure 3) keeps urine bags that are attached to a bed rail firmly in place.


Our **straight IV Pole** can be inserted into the nursing bed patient lifting pole using a plastic adapter. The stand's height is adjustable, and it can support up to 2 kilograms (figure 4).

Our 45-degree **goose-necked IV Pole** prevents wall damage in the event a resident moves laterally while in the Trendelenburg position (large picture).



IV Poles







Additional convenience accessories.

We are committed to making care as comfortable and convenient as possible for residents and their caregivers. To this end, we offer a series of innovative solutions.

The Stiegmeyer bed transporter. This apparatus allows for the transport of nursing beds with fixed wooden frames, and can also be used for bed repair and maintenance purposes. The transporter's large casters swivel individually and have brakes (large picture).

For beds with high mattresses, lateral rail height can be increased by approximately 145 mm using a **bed rail unit** (figure 1) whose tappet is compatible with various rail configurations.

Our **bed rail model 6638** can be mounted on nursing beds that have no track assembly. To mount the product, a threaded sleeve must first be existing on the head and foot elements (figure 2).

Our rotatable and swivellable transmitter holder (figure 3), which features a flexible support, is a highly convenient addition to any nursing bed. The manual switch is simply inserted into the holder (which is not compatible with all of our beds).

The transmitter holder clip (figure 4) can be attached to the upper bed rail, thus allowing the holder to be moved to the most comfortable manual-switch position for the patient (for rails with grasping bars).

Additional convenience accessories.

The **servicing tray** is simply placed on the raised rail, where the tray is very stable. An impact absorbing element under the tray allows for lateral stabilisation. When the tray is not in use, it can be suspended from the headboard (large picture). An armchair model is also available.

The **magazine rack** can be used to store resident records, as well as periodicals, and can be mounted on a bed or bedside table (figure 1).

Our aluminium self-adhesive **name tag holders** allow resident name tags to be removed for marking and then re-inserted (figure 2).

An anti-bedsore device compressor can be suspended from the **device rack**, which is attached to the headboard. The rack can also be used for blankets and towels (figure 3).

The **writing table** is a practical fold-away unit that can be used anywhere and whose position is easy to adjust. The product features a document retention bar at the top edge and aluminium bars at the bottom and sides to prevent documents from sliding off.

The **walker support** can be installed anywhere. Its tappet clamp allows the support to be readily mounted and removed (figure 5).



Serving tray



Bedside table rail





Bedside table accessories.

A bedside table is probably the most important piece of furniture in a resident's room, apart from their bed. We offer a number of accessories that enhance bedside table functionality.

The aluminium **rail** on the top of the bedside table can be used to attach the practical **towel hook** and **bottle holder**, either of which can be slid to any position along the length of the rail (large picture).

The **bottle holder** holds bottles and glasses securely in place on the bedside table (picture 2).

The **bedside table walker holder** is well within the resident's reach, thus making this an optimal mounting position. The holder is firmly attached to the smooth bedside table surface via a tappet clamp (figure 1). The product comes with a self-adhesive mounting plate.

The **power strip** with child-safety element and power cord is VDE-approved and is rated for 230 volt devices (figure 3).

A **covered waste receptacle** with approximately five litres capacity allows waste to be discarded in a safe and hygienic manner.

This plastic receptacle, which is easy to clean, can be mounted on the bedside table's utility bar.

The **magazine** rack can be used to store residents' reading material and frees up space on the bedside table (figure 5).

Stieglmeyer Pflegemöbel GmbH & Co. KG
Pivitsheider Straße 270, 32791 Lage, Germany
Phone: +49 (0) 52 32/98 41 4-30
Fax: +49 (0) 52 32/98 41 4-31
E-mail: care@stieglmeyer.com
www.stieglmeyer.com

EVOCARE AUSTRALIA PTY LIMITED

A.B.N. 98 078 566 604

Trading as **EVOCARE** and **L&M EQUIPMENT**

P.O. Box 1144, Stafford Qld. 4053

Ph: 07 3355 8000 Fax: 07 3355 5043

Website: <http://www.evocare.com.au>

Email: sales@evocare.com.au

workshop@evocare.com.au

warehouse@evocare.com.au

accounts@evocare.com.au



STIEGELMEYER