




# NLMS SPARK SAFE LABORATORY REFRIGERATOR

**NULINE**  
REFRIGERATION

## FEATURES

- Top Mounted System
- E-glass and heated cabinet facia to prevent condensation up to 80% R.H.
- Glass or Solid Door
- Auto Water Disposal
- Hot Gas Defrost
- Castors
- Door Lock
- Digital Display with Built-in Alarms
- Hi/Lo Readout



	2°C to 8°C
	400 L 700 L 1000 L 1614 L
	600w x 620d x 1915h 870w x 700d x 2100h 1200w x 700d x 2100h 1880w x 700d x 2100h



+61 3 9562 3355  
info@nuline.au  
7 Longford Court, Springvale VIC 3171  
[www.nulinerefrigeration.com.au](http://www.nulinerefrigeration.com.au)



# NLMS SPARK SAFE LABORATORY REFRIGERATOR



MODEL	NLMS400/1	NLMS700/1	NLMS1000/2	NLMS1614/3
Gross Capacity	400 L	700 L	1000 L	1614 L
External Dimensions in mm (w/d/h including castors or legs)	600 x 620 x 1915	870 x 700 x 2100	1200 x 700 x 2100	1880 x 700 x 2100
Internal Dimensions in mm (w/d/h)	500 x 470 x 1435	770 x 550 x 1555	1100 x 550 x 1555	1780 x 550 x 1520
Packed Dimensions in mm (w/d/h including pallet)	700 x 720 x 2065	970 x 800 x 2250	1300 x 800 x 2250	1980 x 800 x 2250
Gross Weight – kg (packed weight)	102	112	137	187
Net Weight – kg	95	95	125	175
Voltage	240V	240V	240V	240V
Temperature Range	2°C to 8°C	2°C to 8°C	2°C to 8°C	2°C to 8°C
Ambient Temperature Range	10°C to 32°C	10°C to 32°C	10°C to 32°C	10°C to 32°C
Refrigerant	R134a	R134a	R134a	R134a
Cooling System / Defrost	Static	Static	Static	Static
Exterior Cabinet Finish	White Colourbond	White Colourbond	White Colourbond	White Colourbond
Inner Liner Material	White Colourbond	White Colourbond	White Colourbond	White Colourbond
Door	Solid / Glass E-Glass	Solid / Glass E-Glass	Solid / Glass E-Glass	Solid / Glass E-Glass
Door Hinge	RH	RH	LH / RH	LH / RH / RH
Door Handle	Lockable	Lockable	Lockable	Lockable
Type of Controller	Digital	Digital	Digital	Digital
Temperature Display	External	External	External	External
Data Logger	Optional	Optional	Optional	Optional
Door Ajar Alarms	Optional	Optional	Optional	Optional
High and Low Readout	Visual / Audible	Visual / Audible	Visual / Audible	Visual / Audible
Remote Alarm Contacts (BMS Connectivity for Alarms only) (BMS Connectivity for Temperatures)	Optional Optional	Optional Optional	Optional Optional	Optional Optional
Port Entry	Optional	Optional	Optional	Optional
Battery Back-Up for Controller	Optional	Optional	Optional	Optional
Power Requirement	10amp GPO	10amp GPO	10amp GPO	10amp GPO
Castors	Yes	Yes	Yes	Yes
Adjustable Shelves	5	5	5	15
Usable Shelf Area in mm (w/d)	480 x 450	750 x 530	1080 x 530	565 x 530

Nuline Refrigeration reserves the right to alter the specifications without prior notice.

# NLMS SPARK SAFE LABORATORY REFRIGERATOR



Since 1985, Nuline Refrigeration has been manufacturing refrigerators and freezer cabinets for the Australian market. Our specialised range of biomedical cabinets are used extensively across the hospitals, medical centres, pharmacies and laboratories.

Focusing on customer needs, we manufacture, import and customise the biomedical cabinets from a range of selected global suppliers. Our factory is based in Melbourne, Australia and we are proudly 100% Australian owned business.

Each of our manufactured cabinets are custom built based on customer orders by our craftsmen right here in Melbourne.

We have extensive capability to design, engineer and build a range of refrigerated cabinets for you. Our standard product range includes, fridges, freezers (up to  $-40^{\circ}\text{C}$ ), blanket warmer, mortuary refrigerators and refrigerated incubators.

Distributor

## **EVOCARE AUSTRALIA PTY LIMITED**

A.B.N. 98 078 566 604

Trading as **EVOCARE** and **L&M EQUIPMENT**

P.O. Box 145, Everton Park Qld. 4053

Ph: 07 3355 8000

Website: <http://www.evocare.com.au>

Email: [sales@evocare.com.au](mailto:sales@evocare.com.au)

[workshop@evocare.com.au](mailto:workshop@evocare.com.au)

[warehouse@evocare.com.au](mailto:warehouse@evocare.com.au)

[accounts@evocare.com.au](mailto:accounts@evocare.com.au)



Global-Mark.com.au®

