

# Etac Sling

## Sling with handles

Etac Sling with handles acts as extension of the arms, allowing the caregiver to offer support and assistance whilst maintaining a good ergonomic position and grip, without encroaching on the service user's personal space.

Etac Sling is a flexible tool that allows the caregiver to maintain a more ergonomic grip. It can be used to support many different transfers such as, moving higher up in bed, transfer between sitting positions, etc. Sling has 4 handles on each side. The inside has a non-slip material to help prevent sliding out of position during the transfer, whilst the outside has a low friction surface to aid placement and removal

*Etac Sling is a flexible tool offering a number of different functions.*



**Intended Use**

- Support moving higher up in bed
- Support transfer between bed, wheelchair, chair, etc
- Support stand up / sit down
- Moving back / forward in chair / wheelchair



Description	Width	Length	Item no.
Sling Long, 2 x 4 handles	16 cm 6.3"	190 cm 74.8"	IM428 ✓

**Specifications**

Backside: Nylon  
Padding: Polyester  
Inside: Nylon with polyurethan  
Handles: Polyester



**Washing Instructions**

Do not use fabric softener or bleach.

