






Heavy Duty Laundry Bags

Soiled Linen Bags

Colour Code	Colour Code	Normal Area of Use
R	 RED	Contaminated
Y	 YELLOW	Infectious and Foul
G	 GREEN	Theatre (Sterile Linen)
B	 BLUE	Maternity/Special Items
P	 PURPLE	Chemotherapy
W	WHITE	Ward (General Linen)

**Please specify Bag Colour when ordering

Bag Colours



4H230B-BL, Black

4H230B-B, Blue

4H230B-BR, Brown

4H230B-G, Green



4H230B-GR, Grey

4H230B-P, Purple

4H230B-O, Orange

4H230B-R, Red



4H230B-W, White

4H230B-Y, Yellow

Sharkskin Laundry Bag Features:

- 15kg capacity with lift handles
- 100% heavy duty nylon, cylindrical in shape
- 74cm high x 71cm wide
- Drawstring with non-slip toggle
- Fully reversible
- Pocket for identification

Evocare Soiled Linen Bags are made to Australian Standards and conform to Workplace Health & Safety requirements



Washing Instructions for Laundry Bags:

- Use the same wash cycle used for sheets and a drying time of 4 minutes.
- Laundry Bags are not to be dried with towels as the drying cycle for towels is too long and it will shrink the cord on the bag.