

seca 954

Wireless chair scales to weigh seated patients

EMR integrated



- Patients easy to move from bed to chair due to swiveling arm and foot rests.
- BMI function.
- Smooth running wheels.
- At the press of a button measured results can be transmitted wirelessly to a seca wireless printer or a PC.

seca 954:

Consideration, comfort and safety.

How can one weigh bedridden patients or people with an impaired sense of balance as gently as possible? The answer is a seca chair scale as this can be positioned directly next to the patient, saving the patient – and the carer – unnecessary effort.

Swiveling arm and foot rests for easy transfer to scale.

Consideration and comfort – these are particularly prominent aspects of this stable model with its circular steel tube frame. The armrests which can be folded back make the unimpeded transfer of the patient from chair or bed to the chair scale possible. The foot rests can be folded upwards and pushed under the seat. The lockable rear wheels ensure that it is safe to sit down and stand up.

Effortless mobility.

The seca 954 has four smooth running pivoting wheels making use possible even in the narrowest of spaces.

Clearly laid out control element.

Positioned on the back of the chair scale seca 954 and thus easy to operate is the display and control element with its large LCD display. All clearly laid out functions are easy to control with four large keys. Practical: the carrying case seca 471 provided for storing the power adapter directly on the scale.

Precise results for all cases.

As the zero point of the weighing function can change due to external influences – especially in the case of such a highly sensitive scale, it is automatically checked and aligned electronically once every second. And to ensure that the seca 954 determines the precise weight of restless patients, it is fitted with an effective damping system. An additional function: the BMI value and therefore the nutritional condition of the person weighed can be ascertained automatically at the push of a button. This mobile chair scale, which can be operated by rechargeable battery or power adapter, provides enormous help in rehab, orthopedics, dialysis, nursing homes and pediatrics, allowing small children who cannot yet stand on their own to be weighed while seated.



Paperless and error-free wireless transmission of measurements.

With the seca 360° wireless technology the seca 954 can transmit measurements to the optional seca 360° wireless printer. With the network-capable software solutions seca analytics 105 or seca emr flash 101 and wireless USB adapter seca 456, your PC can receive and analyze measurements and forward them to an Electronic Medical Record (EMR) system. Thus the seca 954 is EMR integrated and ready to handle digital patient records and all the demands the future will bring. More information about EMR integration is at www.seca.com



To make sitting down and standing up easier, the arm and footrests can be folded away.



The ergonomically shaped handles let you push or pull the chair scale in any direction.

seca 954

Technical Data

- Capacity: 300 kg / 600 lbs
- Graduation: 50 g < 150 kg > 100 g / 0.1 lbs < 330 lbs > 0.2 lbs
- Dimensions (WxHxD): 563 x 906 x 970 mm / 22.2 x 35.7 x 38.2"
- Weight: 24.5 kg / 54 lbs
- Power supply: power adapter / rechargeable batteries
- Data transmission: seca 360° wireless technology
- Functions: TARE, pre-TARE, HOLD, auto-HOLD, BMI, auto-CLEAR, kg/lbs switch-over, damping, automatic switch-off, SEND/PRINT
- System compatible with: 360° wireless printer seca 466 advanced, 360° wireless printer seca 465, software seca analytics 105 and seca emr flash 101, 360° wireless USB adapter seca 456

EVOCARE AUSTRALIA PTY LIMITED

A.B.N. 98 078 566 604

Trading as **EVOCARE** and **L&M EQUIPMENT**

P.O. Box 1144, Stafford Qld. 4053

Ph: 07 3355 8000 Fax: 07 3355 5043

Website: <http://www.evocare.com.au>

Email: sales@evocare.com.au

workshop@evocare.com.au

warehouse@evocare.com.au

accounts@evocare.com.au

